



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

## LF00958 - Pediatric Arm

Order code: **4107.LF00958**



Cena bez DPH	714,00 Eur
Price with VAT	863,94 Eur

### Parameters

Cannula, injection, puncturing - filter	Horní končetiny
Quantitative unit	ks

- Practice Venipuncture and Intramuscular Injection Techniques and Procedures Used in Young Children
- Exact Reproduction of the Arm of a Six-Year-Old Child

A soft foam is used to simulate the deltoid muscle and allows the student to "get the feel" of administering intramuscular injections to children. Simulated bone in the shoulder defines and limits the injection area. Water maybe used as an injection fluid in the shoulder area. Includes a Life/form® Arm with white skin, fluid supply bag, 3 cc syringe, one pint Life/form® blood, teaching guide, and hard carrying case. Three-year warranty. Carton size: 22" x 11" x 6" (56 x 28 x 15 cm). Sh. wt. 10 lbs. (4.54 kg).

### Available Supplies

- LF00986 - Replacement Skin and Vein Kits
- LF01056 - Vein Replacements Only
- LF00845 - Venous Blood - Quart. Sh. wt. 0.38 lbs. (0.17 kg)
- LF00846 - Venous Blood - Gallon. Sh. wt. 0.88 lbs. (0.40 kg)