

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz MS-70 - Moisture Analyzer
Order code: 0105.MS70



Cena bez DPH Price with VAT 6.280,00 Eur 7.598,80 Eur

Parameters

Type Balances

Series

Approval (Verif. ČMI)

Capacity (maximum)

Moisture Analyzers

MF, ML, MS, MX

NO

71 g

Min. weighing value0,1Division [%]0,001InterfaceRS-232CCalibrationExternalWeighing pan [mm]Ø 85Quantitative unitks

- Fast and uniform heating with a straight halogen lamp and innovative SRA (Secondary Radiation Assist) technology shorten measurement time
- Standard WinCT-Moisture (for MS-70 & MX-50 only) for easy analyzing mode setting through PC and for real-time graph displaying of moisture content rate change while measuring
- Sodium Tartrate Dehydrate with its stable moisture rate comes standard for accuracy checking
- Optional temperature calibrator, calibration result can be output
- Memory function allows 20 (10 for MF-50, 5 for ML-50) suitable measurement conditions according to sample and 100 (50 for MF-50, 30 for ML-50) measurement result data storage and outputs
- Four choices of measurement programs, Standard, Automatic Complete, Timer, and Manual Mode are selectable
- Clear and easy-to-see, large VFD display shows measurement, setting value, change in moisture content, action status, data number and other important information
- Low maintenance cost assured, thanks to user replaceable halogen lamp without unit downtime with protective chamber for easy clean
- Progress window for heating check
- "Quick Mode" in the measurement program for drying the sample faster by heating the sample at 200°C for 3 minutes.
- Self check function
- Conformity to GLP,GMP,GCP and ISO with date /time ,ID calibration data and check record outputs