

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154

GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

LF04004 - Geri Manikin Auscultation Update

Order code: 4107.LF04004



Information about product price on demand

Parameters

Nursing simulators - filter

Dospělý

Quantitative unit

ks

The GERi™ and KERi™ Auscultation Manikins have all the same features as the Complete GERi™ and KERi™ Manikins with the addition of auscultation sites with heart and lung sounds. For auscultation training, the instructor selects from a menu of 12 heart and lung conditions by a wireless remote control. The student must palpate to identify correct auscultation locations, and will hear different heart and lung sounds as the SmartScope™ is moved to different locations on the manikin. Lung sounds can be detected at five anterior and 12 posterior sites and six anterior heart sites. A diagnosis of condition selected by the instructor can be made by comparing the variations in sounds occurring at the different sites — just like a real person! Each trainer includes manikin, one remote control with LCD display, one SmartScope™ with single-and dual-user headpieces, batteries and hard storage case. Requires two "AA" and two "AAA" batteries (included).

Turn your existing GERi $^{\text{TM}}$ manikin into an Auscultation Trainer. Update package includes torso with auscultation sites, remote control with LCD display, and SmartScope $^{\text{TM}}$.