

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Sakamoto Suction Simulator

Order code: 4109.M175



Cena bez DPH 2.276,00 Eur

Price with VAT 2.753,96 Eur

Parameters

Quantitative unit ks

Basis of respiratory management Endotracheal intubation can be learned.

This model is a simulator with a lifelike texture, and various suction methods can be learned.

The suction of secretions not only in the oral/nasal cavities but also in tubes inserted via the oral cavities is possible.

Since the trachea has already been incised, suction with tracheal cannula insertion can also be performed.

Features

The head can be moved to the left and right.

Practice

Various types of suction training are possible

Suction of secretions in the oral cavity
Suction of secretions in inserted tubes
Suction of secretions in the nasal cavity
Suction of secretions in the tracheal cannula
Washing of the bronchi
Washing of the oral cavity