

## MW1 Pitting Edema Model

Order code: **4103.MW1**



Information about product price on demand

Parameters

Quantitative unit

ks

Five variations of pitting edema models including a normal one. Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th ed.)". The models can be worn on the lower thighs of a SP or a training manikin.

Training items:

- Differentiation of pitting edemas with four different grades and a normal
- Learn the examination procedures for edemas on the lower thighs

Features:

- Four grades of pitting edemas with different resistances, depth of indentations
- Grades of the edemas are based on "Mosby's Guide to Physical Examination (7th Ed.)"
- The models can be worn on the lower thighs of a SP or a training manikin
- The removable fibula plate facilitates training without SP leg or simulation on the soft bases

Degrees of pitting edemas:

- Normal - No indentation left
- Light pitting (+) - Indentations with approx. 2 mm depth
- Moderate pitting (2+) - Indentations with approx. 4 mm depth
- Deep pitting (3+) - Veins, bones and joints are indistinct because of swell, indentations with approx. 6 mm depth
- Very deep pitting (4+) - Indentations with approx. 8 mm depth

Set includes:

- Strap-on pitting edema model 5 variations (1 piece each)