



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

LIFEPAK 1000 with ECG display

Order code: **4401.99425000126**



Information about product price on demand

Parameters

Quantitative unit

ks

The LIFEPAK 1000 defibrillator is the most rugged AED. ADAPTIV biphasic technology provides a range of energy up to 360 joules, and an upgradeable platform. The 1000 is equipped with cprMAX™ technology, enabling care providers to change protocols as standards of care evolve. A large, intuitive screen displays graphics and ECG readings that are clear and easy to read from a distance. A visible and accessible battery gauge gives rescuers the ability to assess battery readiness. The non-rechargeable battery can deliver more than 425 shocks.

LIFEPAK 1000 defibrillator

You can focus on the patient, on the moment, on what matters.

The LIFEPAK 1000 defibrillator ensures you have the best equipment to fulfill the commitment you made when you decided to enter this field—to save lives.

Although not everyone can be saved from sudden cardiac arrest, studies show that early defibrillation can dramatically improve survival rates.

Customized Care for Your Patient

- Easy customization of settings to meet your protocol needs in the field
- cprMAX™ technology feature gives you flexibility, so you can increase hands-on time with your patient
- ADAPTIV™ biphasic technology provides a full range of energy from 150 to 360 joules

A Rugged and Intuitive Defibrillator

- The most rugged AED ever built by Physio-Control, with an IP55 rating
- Visible and accessible battery gauge gives you the ability to assess battery readiness
- Battery performs advanced communications directly with the device so you always know the battery's condition
- Large intuitive screen displays graphics and ECG readings Manual defibrillation (AED override) option is available

A Resuscitation Platform for Now and the Future

- Upgradeable platform and flexible setup options to help you focus on continually improving patient outcomes
- The 1000 is compatible with other Physio-Control products and services across the full spectrum of care, from the field to the hospital, including data download

AED users should be trained in CPR and the use of the AED. Please consult your physician. LIFEPAK AEDs require a prescription.