

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Tooth Extraction Model SUG2004-UL-SP-

DM-28

Order code: 4112.342577



Cena bez DPH 237,00 Eur
Price with VAT 286,77 Eur

Parameters

Quantitative unit ks

- Tooth extraction practice can be done on 6 maxillary teeth and 6 mandibular teeth
- Anatomically shaped extraction teeth are held in place in the socket by hot-melt (heated glue) with the right amount of resistance allowing realistic extraction
- Extraction pressure, direction and feel can be experienced at a more clinically accurate level
- A more advanced practice can be done by attaching the model to a phantom or mannequin
- Repeated practice can be done by using a special hot-melt glue gun (*please inquire for more information)
- This model cannot be used with the FE articulator (prices are for jaws with a D articulator)