



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

B-bag - neonatal
Order code: **4410.BG07021**



Information about product price on demand

Parameters

Accessories and aids

Vaky

The double wall structure allows resuscitation professionals to monitor the airway pressure of the patient with the fingertips, revealing with timeliness any possible pressure peaks which are limited by the double wall. All of the connections conform to the specifications of the ISO. The resuscitator is light enough to be sustained and used for long periods of time without fatigue. In cases of necessity the patient valve can be connected to the resuscitator by a corrugated tube. Available in the following versions:

- **Neonatal:** for patients below 7 kg

Single packaging in cardboard box printed using a four-coloured process.