

**AirSim Combo Bronchi X**  
Order code: **4110.CTC95100X**



Information about product price on demand

Parameters

Quantitative unit

ks

TruCorp AirSim Combo Bronchi X features the uniquely constructed AirSim® X airway with enhanced with palpable cricoid landmarks, laryngeal cartilages and tracheal rings. Perfect for practicing intubation, ventilation, percutaneous tracheostomy, needle and surgical cricothyroidotomy plus diagnostic bronchoscopy techniques including lung suctioning & isolation.

### **Model Features**

- Innovative AirSim® X airway and nasal cavity created from CT DICOM data
- Anatomically correct internal features and visually accurate landmarks in airway & nasal passage
- Wraparound replaceable neck skin rotates for 10-15 incisions
- Easy identification of sternal notch, tracheal rings & cricoid & laryngeal cartilages
- Inflatable tongue with lifelike texture can simulate edema
- Silicone skin covering with lifelike texture
- Replaceable lung bags inflate with successful ventilation

### **Medical Procedure Training**

- Needle and surgical cricothyroidotomy
- Percutaneous tracheostomy
- Diagnostic bronchoscopy
- Lung suctioning techniques
- Lung isolation techniques using left & right endobronchial tubes & bronchial blockers
- Double nasotracheal intubation
- Bag-valve mask (BVM) ventilation techniques
- Full range of supraglottic device insertion
- Direct and video laryngoscopy including Sellick Maneuver
- Awake fiberoptic examination
- Endotracheal tube insertion
- Combi tube insertion
- Naso gastric tube insertion techniques

### **Package Contents**

- 1 Carrier/Storage Case
- 1 Bottle Of TruCorp® Lubrication
- 3 Adult Lung Bags replicating 2 lungs and stomach
- 5 Larynx Inserts (1 attached, 4 extra)
- 1 Wrap Around Neck Skin (attached)
- 2 sets of adhesive pads (attachable to suction cups for increased stability if required)

### **Parts & Accessories**

- AirSim Combo X Replacement Airway
- Replacement Head Skin And Nasal Passage
- Base Plinth
- Neck Locks
- Suction Cup (Set Of 8)