


## 8000973 - Medtronic Capnostream™ 35 Patient Monitor Screen Simulation for REALITi360

Order code: **4127.8000973**



Information about product price on demand

Parameters

Quantitative unit

ks

The simulated patient monitor “Capnostream™ 35” is an optional premium branded screen simulation that works with the REALiTi360 vital sign monitor and defibrillator app. It is designed to be used with the patient monitor simulation systems REALiTi Plus and REALiTi Pro by iSimulate.

Instructors can easily add training scenarios or download new scenarios from the REALiTi360 COMMUNITi platform.

REALiTi360 Premium Screens mimic the user interface of several leading patient monitor and defibrillator manufacturers. Through collaboration with popular defibrillator and monitor manufacturers, iSimulate has created hyper-realistic interfaces which look and work like the real vital sign monitors and defibrillators used at the point of care. The current and growing selection of patient monitors and defibrillators allows students to experience a significant number of monitors and defibrillators that they might encounter in their EMT or Hospital careers.

Based on individual training needs and the product configuration, users can select one or more of these Premium Display Screens for their simulation trainings.

The patient monitor simulation of “Capnostream™ 35” gives REALiTi360 users the option to train respiratory monitoring of patients in any clinical setting. This screen simulation perfectly mimicks the “Capnostream™ 35” monitor including all of its branded features.

### **Required**

- 4127.1022815 - Simulated Patient Monitor - REALiTi Plus
- 4127.1022816 - Simulated Patient Monitor with Debriefing & CPR Feedback - REALiTi Pro